


# Operating Table

## DST-500A

Great performance in a compact size

OPERATING ROOM



Modern OR room amazing Operating Table--- DST-500A. Application of modern technology, combined with the idea of humanized design, meet the requirements of operation in the process of video observation and diagnosis, comprehensive strength, for each part of the human body surgery hospital operating room.

DST-500A --It bears patient weight up to 280kg ( 617 lbs), with a super lower position comes to 500mm(19.69in), permits ergonomic working positions, both seated and standing.

Except as an universal operating table, DST-500A has the function X-ray imaging, match up to C-arm, which is economic to modern demands.

# Operating Table

## DST-500A

Great performance in a compact size

OPERATING ROOM

### Technical Data

Bed Panel Size	2040 x 520mm
Bed height from the ground (without pads)	Min height: 850 mm Max height: 1050mm
Lifting route (electric)	200 ± 20mm
Horizontal route (electric)	300 ± 20mm
head and foot tilt (electric)	head tilt: 20 ° ± 2 °
	foot tilt: 20 ° ± 2 °
Lateral tilt (electric)	Left: 19° ± 2 °
	Right: 19° ± 2 °
leg plate fold (manual)	downward: 90 °
	Horizontal extend: 90 °
back plate fold (electric)	upward: 65°
	Downward: 30°
Head section fold (manual)	Upward: 20° ± 2 °
	Downward: 90° ± 2 °
External kidney bridge lifting route(manual)	120 ± 10mm
Net weight/ capacity weight	170kg/ 300kg
Table locking	manual
Power supply	220V 50Hz 750VA

### Features

1. It is an economic one for universal operating table, various surgeries.
2. Table board: dedicated composite material board with high X-ray tranmittancy and favorable uniformity, satisfy the imaging demand.
3. Built in kidney bridge, match up to kidney, cervical and lumbar surgery.
4. Stable and noiseless electric control for lifting up and down, trendelenburg and anti-trendelenburg, turning right and left.
5. Seperated leg board: The leg board can be disassembled and folded, apply to gynecology and urology surgeries.
6. Central braking casters: moveable and steady.

### Standard Accessories



Y-001  
mattress pad




Y-002  
cable controller



Y-003  
lateral support



P-001/P-002  
lock block/clamping block



Y-004  
carbon fiber arm support



Y-006 fixed belt



Y-009 leg support



Y-008 anaesthesia screen



D-001 orthopedic tractor