

HYDRAULIC PET OPERATING TABLE

SUS304# stainless steel table top

GST-100 series



Model: GST-100


Hydraulic Veterinary Operating Table



Decline 65cm

Pedal Hydraulic

Pedal Hydraulic Lifting Without Power
Convenient Work Anytime And Anywhere.



304# stainless steel
countertops

Iv stand

Appliance plate

Liquid waste barrels

Foot hydraulic

Fixed foot

304 # Stainless Steel Panel

*Thickening stainless steel panel
surface is smooth
easy to clean
thicken panel
bearing ability*



304 # stainless steel Bracket

*Adopts high quality stainless steel
with heat resistant
high temperature resistant*



Waste Barrels/ Sewage Discharge

Waste barrels of optional disassembly
convenient cleaning/sewage discharge
convenient and practical



IV Stand/Tray

height size can adjust at will
tray removable,
stainless steel materials



size can be customized

Rise to the top

108cm

Drop to the bottom

30cm

Supply Two Kinds Of table Wheels

Universal wheel



Adjusting wheel



&

Product Description

Pet lifting operating table, hydraulic, electric, pedal, constant temperature, medical 304 stainless steel diagnosis and treatment table, non-slip, firm, stable, antibacterial, suitable for pet treatment, sterilization, caesarean section, rehabilitation physiotherapy, pet grooming, etc.

specification

table length	120-132mm
table width	600-670mm
height adjust	30-108cm
material	304medical stainless steel
lifting method	Hydraulic/ electric
size	small and large
small size	120*60*65*95cm
small size packing	124*63*73cm
large size	132*67*65*100cm
large size packing	132*71*73cm